



TURQUOISE BIOMEDICAL TECHNOLOGIES

*Manufactures
Happiness! :)*

PLUSHINE

System Accessories for Optimizing of Details

www.tbmedical.com


+ Solutions for an efficient workflow

PLUSHINE series system accessories are at your service for the entire flow of goods in a CSSD. Designed for your different needs for pre-cleaning, inspection, packing, storage and distribution; PLUSHINE series system accessories are ideal partner for increased efficiency. You can always combine your PLUSHINE series system accessories with your PLUSEAL series sealers or with PLUSONIC series of ultrasonic cleaners for increased efficiency in workflow. Even though they are carefully designed to be flexible, we are ready to serve you for your specific needs for different designs and sizes.

+ Wide range of selection

PLUSHINE series system accessories combine your needs from wide range of containers & baskets to storage & goods carrying solutions and even further to support devices such as spray gun rinser and manual trolley washing station.

Manual trolley washing station is a time saving solution for achieving hygiene in the trolley used in hospital. The flexibility of system is increased with 15 m. long hose, combined with hose hanger. 10 L/min. flow rate for water, which is diluted with chemical from one of the 15 adjustable dilution ratios offer a satisfying solution. Manual trolley washing station can operate up to 8,6 bar and 71 °C, thanks to its robust stainless steel design.



+ PLUSHINE series stainless equipments with wide range of standardized solutions and optional features for tailored solutions will make you feel the comfort in process. AISI 304 grade stainless steel, with perfect finishing, increases the robustness of your selection of stainless equipments. To offer you case specific solutions, we can assist you for your custom needs. Please see our system accessories list below and don't hesitate to contact us for detailed information



Please see our system accessories list below and don't hesitate to contact us for detailed information

PLUSHINE SERIES SUPPORTING DEVICES

ORDER NUMBER	DESCRIPTION	DIMENSIONS (W x D x H)
711001	MANUAL TROLLEY WASHING STATION	203 x 137 x 235 mm.
711002	SPRAY GUN RINSER W/ 8 NOZZLES	-

PLUSHINE SERIES STAINLESS EQUIPMENTS

ORDER NUMBER	DESCRIPTION	DIMENSIONS (W x D x H)
713001	WORKSTATION FOR C.S.S.D (PACKING TABLE)	1150 x 800 x 900 mm.
711002	ROLL HOLDER	858 x 425 x 200 mm.
711003	ROLL CUTTER	858 x 425 x 200 mm.
711004	WORKTABLE WITH STANDART WORKSURFACE ONLY	900 x 650 x 900 mm.
711005	WORKTABLE WITH STANDART WORKSURFACE ONLY	1150 x 650 x 900 mm.
711006	WORKTABLE WITH STANDART WORKSURFACE ONLY	1400 x 650 x 900 mm.
711007	WORKTABLE WITH STANDART WORKSURFACE ONLY	1900 x 650 x 900 mm.
711008	WORKTABLE WITH STANDART WORKSURFACE ONLY	2400 x 650 x 900 mm.
711009	WORKTABLE WITH STANDART WORKSURFACE ONLY	2800 x 650 x 900 mm.
711010	WORKTABLE WITH BOTTOM SHELF	900 x 650 x 900 mm.
711011	WORKTABLE WITH BOTTOM SHELF	1150 x 650 x 900 mm.
711012	WORKTABLE WITH BOTTOM SHELF	1400 x 650 x 900 mm.
711013	WORKTABLE WITH BOTTOM SHELF	1900 x 650 x 900 mm.
711014	WORKTABLE WITH BOTTOM SHELF	2400 x 650 x 900 mm.
711015	WORKTABLE WITH BOTTOM SHELF	2800 x 650 x 900 mm.
711016	WORKTABLE WITH BOTTOM & MIDDLE SHELVES	900 x 650 x 900 mm.
711017	WORKTABLE WITH BOTTOM & MIDDLE SHELVES	1150 x 650 x 900 mm.
711018	WORKTABLE WITH BOTTOM & MIDDLE SHELVES	1400 x 650 x 900 mm.
711019	WORKTABLE WITH BOTTOM & MIDDLE SHELVES	1900 x 650 x 900 mm.
711020	WORKTABLE WITH BOTTOM & MIDDLE SHELVES	2400 x 650 x 900 mm.
711021	WORKTABLE WITH BOTTOM & MIDDLE SHELVES	2800 x 650 x 900 mm.
711022	WATER FILTRATION TABLE	600 x 600 x 900 mm.
711023	INSTRUMENT TRANSPORT TROLLEY	450 x 450 x 840 mm.
711024	INSTRUMENT TRANSPORT TROLLEY	610 x 450 x 840 mm.
711025	INSTRUMENT TRANSPORT TROLLEY	750 x 450 x 840 mm.
711026	INSTRUMENT TRANSPORT TROLLEY	910 x 450 x 840 mm.
711027	INSTRUMENT TRANSPORT TROLLEY	750 x 610 x 840 mm.
711028	INSTRUMENT TRANSPORT TROLLEY	910 x 750 x 840 mm.
711029	INSTRUMENT TRANSPORT TROLLEY	1200 x 610 x 840 mm.
711030	CONTAINER CARRIAGE TROLLEY	850 x 500 x 1550 mm.
711031	CONTAINER CARRIAGE TROLLEY	1030 x 630 x 1247 mm.
711032	CLOSED CONTAINER CARRIAGE TROLLEY	1150 x 600 x 1621 mm.
711033	CLOSED CONTAINER CARRIAGE TROLLEY	1195 x 1050 x 1621 mm.
711034	PAPER DISPENSING TROLLEY	1350 x 600 x 845 mm.
711035	SLIPPER TROLLEY	1000 x 540 x 1770 mm.
711036	STEP STOOL WITH SINGLE STEP	600 x 350 x 120 mm.
711037	STEP STOOL WITH SINGLE STEP	600 x 350 x 220 mm.
711038	STEP STOOL WITH SINGLE STEP	600 x 350 x 300 mm.
711039	STEP STOOL WITH SINGLE STEP	600 x 650 x 400 mm.
711040	WALL CABINET (MATERIAL CUPBOARD)	450 x 900 x 1800 mm.
711041	WIRE BASKETS TROLLEY	1000 x 540 x 1770 mm.
711042	INSTRUMENT WASHING SINK WITH 1 SINK	1200 x 600 x 900 mm.
711043	INSTRUMENT WASHING SINK WITH 1 SINK	1400 x 600 x 900 mm.
711044	INSTRUMENT WASHING SINK WITH 2 SINKS	1400 x 600 x 900 mm.
711045	INSTRUMENT WASHING SINK WITH 2 SINKS	1600 x 600 x 900 mm.
711046	INSTRUMENT WASHING SINK WITH 2 SINKS	1800 x 600 x 900 mm.
711047	INSTRUMENT WASHING SINK WITH 3 SINKS	1800 x 600 x 900 mm.
711048	INSTRUMENT WASHING SINK WITH 2 SINKS	2000 x 600 x 900 mm.
711049	INSTRUMENT WASHING SINK WITH 3 SINKS	2000 x 600 x 900 mm.

ORDER NUMBER	DESCRIPTION	DIMENSIONS (W x D x H)
713050	INSTRUMENT WASHING SINK WITH 2 SINKS	2400 x 600 x 900 mm.
711051	INSTRUMENT WASHING SINK WITH 3 SINKS	2400 x 600 x 900 mm.
711052	LINEN INSPECTION AND FOLDING TABLE	1400 x 750 x 900 mm.
711053	LINEN INSPECTION AND FOLDING TABLE	1400 x 1400 x 900 mm.
711054	LINEN INSPECTION AND FOLDING TABLE	1800 x 750 x 900 mm.
711055	LINEN INSPECTION AND FOLDING TABLE	1800 x 1400 x 900 mm.
711056	LINEN INSPECTION AND FOLDING TABLE	2400 x 750 x 900 mm.
711057	LINEN INSPECTION AND FOLDING TABLE	2400 x 1400 x 900 mm.
711058	WIRE STORAGE SHELF WITH 4 LEVELS	610 x 360 x 1600 mm.
711059	WIRE STORAGE SHELF WITH 4 LEVELS	1220 x 360 x 1600 mm.
711060	WIRE STORAGE SHELF WITH 4 LEVELS	1830 x 360 x 1600 mm.
711061	WIRE STORAGE SHELF WITH 5 LEVELS	610 x 360 x 1830 mm.
711062	WIRE STORAGE SHELF WITH 5 LEVELS	1220 x 360 x 1830 mm.
711063	WIRE STORAGE SHELF WITH 5 LEVELS	1830 x 360 x 1830 mm.
711064	CONTROLLED PASS-THRU DOOR SYSTEM	-
711065	WALL SHELF WITH SINGLE LEVEL SHELF	900 x 300 x 300 mm.
711066	WALL SHELF WITH DOUBLE LEVEL SHELVES	900 x 300 x 500 mm.
711067	CLOSED CARRIAGE TROLLEY	600 x 400 x 730 mm.
711068	CLOSED CARRIAGE TROLLEY	900 x 640 x 745 mm.
711069	CLOSED CARRIAGE TROLLEY	900 x 640 x 1045 mm.
711070	CLOSED CARRIAGE TROLLEY (6 STU CAPACITY)	995 x 640 x 745 mm.
711071	CLOSED CARRIAGE TROLLEY (9 STU CAPACITY)	995 x 640 x 1045 mm.

PLUSHINE SERIES BASKETS

ORDER NUMBER	DESCRIPTION	DIMENSIONS (W x D x H)
711072	1 STU BASKET	575 x 280 x 265 mm.
711073	1/2 STU BASKET	575 x 280 x 135 mm.
711074-71084	VARIOUS SIZES OF WIRE BASKETS	-
711085	WIRE BASKET - CUSTOM SIZE	-
711086	ENDOSCOPY WIRE BASKET - S SIZE	290 x 80 x 55 mm.
711087	ENDOSCOPY WIRE BASKET - L SIZE	460 x 80 x 55 mm.
711088	STANDARD MESH BASKET - S SIZE	255 x 245 x 50 mm.
711089-711119	VARIOUS SIZES OF MESH BASKETS	-
711120	MESH BASKET - CUSTOM SIZE	-

PLUSHINE SERIES CONTAINERS & ACCESSORIES

ORDER NUMBER	DESCRIPTION	DIMENSIONS (W x D x H)
715001 - 715030	S SIZE CONTAINERS - BOTTOM NON PERFORATED/PERFORATED - STANDART/GUARDING/BARRIER	285 x 280 x - mm.
715031 - 715048	M SIZE CONTAINERS - BOTTOM NON PERFORATED/PERFORATED - STANDART/GUARDING/BARRIER	465 x 280 x - mm.
715049 - 715078	L SIZE CONTAINERS - BOTTOM NON PERFORATED/PERFORATED - STANDART/GUARDING/BARRIER	580 x 280 x - mm.
715079 - 715084	MINI CONTAINERS - BOTTOM NON PERFORATED/PERFORATED - S/M/L SIZE	300 x 140 x - mm.
715085	IMPLANT CONTAINER - BOTTOM PERFORATED - M SIZE	310 x 138 x 65 mm.
715086	IMPLANT CONTAINER - BOTTOM PERFORATED - L SIZE	500 x 155 x 75 mm.
711087	SCOPY CONTAINER - BOTTOM PERFORATED - M SIZE	250 x 60 x 50 mm.
711088	SCOPY CONTAINER - BOTTOM PERFORATED - L SIZE	450 x 70 x 70 mm.
711089-711096	DENTAL CONTAINERS - BOTTOM NON PERFORATED/PERFORATED - S/M/L/XL SIZE	310 x 190 x - mm.
719001 - 719006	MESH BASKET FOR MINI/DENTAL CONTAINERS - S/M/L SIZE	265 x 125 x 30 mm.
719007 - 719027	DIVIDERS FOR S/M/L SIZE CONTAINERS - XS/S/M/L/XL SIZE	-
719028 - 719031	WIRE BASES FOR MINI, S, M & L SIZE CONTAINERS - P/S/M/L SIZE	-
719032 - 719037	SILICONE MATS - XXS/XS/S/M/L/XL SIZE	-
719038 - 719043	LABELS - GRAY/YELLOW/RED/GREEN/BLUE/BLACK	-
719044	PAPER LABEL W/INDICATOR FOR S, M & L CONTAINERS	-
719045	PAPER LABEL W/INDICATOR FOR DENTAL & MINI CONT.	-
719046 - 719057	SECURITY SEALS, CLAMPS, FILTERS & FILTER HOLDERS	-

This document is intended for the use of partners of TURQUOISE BIOMEDICAL TECHNOLOGIES Inc. Reproduction in whole or in part by any other party is prohibited. TBT specifications are given for reference only and are subject to change without any previous notice. Accessory/ Product images are for illustrative purposes only. rev:EN 00001-2013 © 2013, TBT Inc. All rights reserved.



TURKUJAZ BIYOMEDİKAL TEKNOLOJİLER VE SAĞLIK HİZMETLERİ SAN. TİC. LTD. STİ.
 İkitelli Org. San. Sitesi İsteksan Sitesi D1 Blok No:10 İkitelli - Basaksehir / İstanbul
 Tel: +90 212 320 37 53 Fax: +90 212 320 53 51 info@tbtmedical.com - www.tbtmedical.com