

TECHNICAL DATASHEET FOR PLUSTEAM SERIES SILENT AIR COMPRESSOR

Reference

519002 519004

Designation

Silent Air Compressor

Description

Silent air compressor system for internal use in CSSD.



Usage

PLUSTEAM Silent Air System is used to supply compressed air for operating valve system and doors.

Quality Certification & Compliance

All process across **TBT** is run in compliance with EN ISO 9001 & EN ISO 13485.

Silent air compressor is not subject to CE.

Materials Used & Main Features

Stainless covering with painting on is forming the main body. Unit is furnished with 2 wheels in order to enable mobility. A grab handle and 2 standing feet with plastic bottoms are placed on unit.

Available Types & Components

Ref. No	Power	Air Inlet	Pressure	Tank Volume	Sound Level	Weight
519002	700 W (1 HP)	150 L/min	8 bar (115 psi)	24 L	50 dB (A)	21 kg
519004	1000 W (1,5 HP)	196 L/min	8 bar (115 psi)	40 L	62 dB (A)	38 kg

Please request information from your Agent / Sales Representative for getting more information on tailor made solutions of **PLUSTEAM** series.

Cleaning

Clean the unit with a dry fabric. Use only stainless steel maintenance spray for upgraded level of protection.



Comments

The product must be used by specialised personnel only. Check the package's integrity before use. Do not use if the package is open or damaged.

Our company has various designs and several technical variants. There may be differences between the same product types and various versions. This document is intended for the use of Distributors/Dealers/Customers of TBT Inc. Reproduction in whole or in part by any other party is prohibited. Specifications are given for reference only and are subject to change without any previous notice. Accessory/Product images are for illustrative purposes only. rev.EN 003-2017 © 2017, TBT Inc. All rights reserved.