




Manufactures Happiness! :)

PLUSHINE SERIES SYSTEM ACCESSORIES DATAFILE

STAINLESS EQUIPMENTS

MODEL	CAT. NO	DESCRIPTION	DIMENSIONS (W x D x H)
	713003	ROLL CUTTER	858 x 425 x 200 mm.
	<ul style="list-style-type: none"> - For cutting sterilization reels with its moving blade - Bar for holding sterilization reels - Made of AISI 304 stainless steel 		
	<p>OPTIONALS</p> <p>712024 WALL HANGER WITH 2 X SHELVES, 1 X METAL CUP & 1 X DRAWER</p>		